


## AUXILIARY CLASSIC ELECTRODE ACE.E\*



### Description:

ACE is a glass tube with **platinum plate wire**. Standard cable length of auxiliary classic electrode is 1,5 m. Connector or banana plug color is black. Auxiliary classic electrode (ACE) is designed for measurements with thermostated cell TC3.

### Dimensions:

D1	D2	Φ1	Φ2
60 mm	40 mm	10 mm	3 mm

Weight: 24 gms

### Electrode material:

#### ACE.E\* (Auxiliary Classical Electrode)


Electrode is defined when replacing asterisk by letter or number.

E – Electrode ending

S... bare wire

1...banana plug

4... BNC connector



### Ordering Information

- The order is specified by whole product code
- Minimum order quantity – 1 electrode
- Delivery time for standard ACE.E\* electrode is 4 weeks from receipt of order
- Delivery time for non-standard ACE.E\* electrode depends on final technical specification of order

Example of Order:

- 10 pieces – ACE.E1