

## CONNECTOR FOR ELECTROCHEMICAL SENSORS TYPE AC5

Type: KA5.\*


### Description:

The connector enables the use of the electrochemical sensor AC5.

### Physical Parameters:

#### *Dimensions:*

Weight: 100 gms  
 Height: 7 mm  
 Length: 50 mm  
 Width: 30 mm



### Types of termination


Model	Cable	Termination
KA5.S	Three core Shielded cable	Single conductors
KA5.1	Three core Shielded cable	Banana plugs
KA5.2	Three core Shielded cable	BNC connectors

### Example of Order

- 5 pieces - KA5.S

### Evaluating Units

- Any potentiostat with 8 channels (check the cable termination to fit to your potentiostat)



### Ordering Information

- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA5 connector is 4 weeks from receipt of order
- Delivery time for non-standard KA5 connector depends on final technical specification of order