

# THERMOSTATED CELL

Type : TC2



## Description

Simple glass cell for electrochemical measurements jacket which enables the thermostating.

## Physical Parameters

Dimensions		
Length	Height	Width
70 mm	95 mm	70 mm

Weight = 155 gms

## Technical Parameters

Volume A: 10 ml

Volume B: 50 ml

## Optional Equipment

- MFS.\*

## Device Usage

- Electrochemical measurements
- Measurements with biosensors

## Ordering Information

- The order is specified by whole product code
- Minimum order quantity – 1 thermostated cell
- Delivery time for standard TC2 is 4 weeks from receipt of order
- Delivery time for non-standard TC2 depends on final technical specification of order

Example of Order:

- 10 pieces – TC2

