

# THERMOSTATED CELL

Type : TC3



## Description

Simple jacketed glass cell for electrochemical measurements with miniaturised classical electrodes WCE, RCE and ACE. The cell has a lid with four holes to fix electrodes in the cell.

## Physical Parameters

Dimensions		
Length	Height	Width
70 mm	80 mm	70 mm

Weight = 219 gms


## Technical Parameters

Volume A: 10 ml

Volume B: 30 ml

## Optional Equipment

- ACE.\*
- RCE.\*
- WCE.\*



## Device Usage

- Electrochemical measurements

## Ordering Information

- The order is specified by whole product code
- Minimum order quantity – 1 thermostated cell
- Delivery time for standard TC3 is 4 weeks from receipt of order
- Delivery time for non-standard TC3 depends on final technical specification of order

Example of Order:

- 10 pieces – TC3