


**SMALL CONNECTOR FOR ELECTROCHEMICAL SENSORS  
 TYPE AC8**

---

**Type: KA8s.\***

**Description**

The connector enables the use of the biosensor based on the substrates AC8 in classical electrochemical arrangement.



**Physical Parameters**

*Dimensions:*

Weight: 35 gms  
 Length: 6 mm  
 Width: 8 mm

**Types of termination**

Model	Cable	Termination	Evaluating Units
KA8.S	Four core Shielded cable	Single conductors	any device
KA8.1	Four core Shielded cable	Banana plugs	any device
KA8.4	Four core Shielded cable	BNC connectors	any device
KA8.5	Four core Shielded cable	Small banana plugs	any device

**Ordering Information**


- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA8s connector is 4 weeks from receipt of order
- Delivery time for non-standard KA8s connector depends on final technical specification of order

**Example of Order**


- 5 pieces - KA8s.S

**Internal Wiring**


KA8s.S



KA8s.1



KA8s.4



KA8s.5

