

CONNECTOR FOR ELECTROCHEMICAL SENSORS TYPE AC9

Type: KA9.S


Description

The connector enables the use of the electrochemical sensor AC9.

Physical Parameters

Dimensions:

Height: 33 mm
 Length: 36 mm
 Width: 23 mm



Types of termination

Model	Cable	Termination	Evaluating Units
KA9.S	Nine core Shielded cable	Single conductors	Any potentiostat with 8 channels (check the cable termination to fit to your potentiostat

Ordering Information

- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA9 connector is 8 weeks from receipt of order
- Delivery time for non-standard KA9 connector depends on final technical specification of order

Example of Order

- 5 pieces - KA9.S