

AUXILIARY CLASSIC ELECTRODE

Type: ACE.E*

Description

ACE is a glass tube with **platinum plate wire**. Standard cable length of auxiliary classic electrode is 1,5 m. Connector or banana plug color is black. Auxiliary classic electrode (ACE) is designed for measurements with thermostated cell TC2, TC3.



Physical Parameters

Dimensions:

Length: 57 mm
 Diameter: 11 mm
 Weight: 24 gms

Electrode material

ACE.E* (Auxiliary Classical Electrode) - Pt

A - Auxiliary	1 - Banana plug
C - Classical	2 - 7 poles BVT connector
E - Electrode	3 - Triad 01 PalmInstruments
E - Electrode ending	4 - BNC connector
S - Bare wire	5 - Small banana plugs


Ordering Information

- The order is specified by whole product code
- Minimum order quantity - 1 electrode
- Delivery time for standard ACE.E* electrode is 4 weeks from receipt of order
- Delivery time for non-standard ACE.E* electrode depends on final technical specification of order


Example of Order - 10 pieces - ACE.E1

Internal Wiring


ACE.E5



ACE.E1



ACE.E2



ACE.E3

ACE.E4



ACE.E5

