


## WORKING CLASSIC ELECTRODE

Type: WCE.W\*.D\*\*.E\*\*\*

### Description

WCE is a glass tube with a Pt or Au wire inside. Only circle glossy area of the metal wire is outside the glass tube. Standard cable length of working classic electrode is 1,5 m. Connector or banana plug color is red. Working classic electrode (WCE) is designed for measurements with TC2, TC3 glass cell.



### Physical Parameters

#### *Dimensions:*

Length: 57 mm  
 Diameter: 11mm  
 Weight: 50 gms

### Electrode material

WCE.W\*.D\*\*.E\*\*\* (Working Classical Electrode)

W - Working	3 - platinum wire, $\Phi$ 1000 $\mu$ m
C - Classical	4 - golden wire, $\Phi$ 20 $\mu$ m
E - Electrode	E - Electrode ending
W - Electrode material	S - Bare wire
S - Pt	1 - Banana plug
1 - Au	2 - 7 poles BVT connector
D - Diameter of the wire	3 - Triad 01 PalmInstruments
1 - platinum wire, $\Phi$ 100 $\mu$ m	4 - BNC connector
2 - platinum wire, $\Phi$ 600 $\mu$ m	5 - Small banana plugs

### Ordering Information


- The order is specified by whole product code
- Minimum order quantity - 1 electrode
- Delivery time for standard WCE.W\*.D\*\*.E\*\*\* electrode is 4 weeks from receipt of order
- Delivery time for non-standard WCE.W\*.D\*\*.E\*\*\* electrode depends on final technical specification of order

Example of Order - 10 pieces - WCE.WS.D1.E4

**Internal Wiring**


**WCE.W\*.D\*\*.ES**

**WCE.W\*.D\*\*.E1**



**WCE.W\*.D\*\*.E2**

**WCE.W\*.D\*\*.E3**



**WCE.W\*.D\*\*.E4**

**WCE.W\*.D\*\*.E5**

