

AUXILIARY CLASSIC ELECTRODE CONIC

Type: ACEc.E*

Description

ACEc is a glass tube with **platinum plate wire**. Standard cable length of auxiliary classic electrode is 1,5 m. Connector or banana plug color is black. Auxiliary classic electrode (ACEc) is designed for measurements with thermostated cell TC4, TC5, TC6.

Physical Parameters

Dimensions:

Length: 57 mm
 Diameter: 11 mm
 Weight: 24 gms

Electrode material

ACEc.E* (Auxiliary Classical Electrode Conic) - Pt

A - Auxiliary	1 - Banana plug
C - Classical	2 - 7 poles BVT connector
E - Electrode	3 - Triad 01 PalmInstruments
c - Conic	4 - BNC connector
E - Electrode ending	5 - Small banana plugs
S - Bare wire	


Ordering Information

- The order is specified by whole product code
- Minimum order quantity - 1 electrode
- Delivery time for standard ACEc.E* electrode is 4 weeks from receipt of order
- Delivery time for non-standard ACEc.E* electrode depends on final technical specification of order


Example of Order: - 10 pieces - ACEc.E1

Internal Wiring


ACEc.ES



ACEc.E1



ACEc.E2



ACEc.E3



ACEc.E4



ACEc.E5

