

**CONIC CONNECTOR FOR ELECTROCHEMICAL SENSORS
 TYPE AC1, AC4, AC11, AP1, CC1, CC2, TS1**

Type: KA1C.*

Description

The connector enables the use of the biosensor based on the substrates AC1, AC4, AC11, AC13, AP1, CC1, CC2 and TS1 in classical electrochemical arrangement with stirred vessel. The cell TC4, TC5, TC6 is recommended for use with KA1C.*.



Physical parameters

Dimensions:

Weight: 36 gms
 Length: 37 mm
 Width: 10 mm

Types of termination

Model	Cable	Termination	Evaluating Units
KA1C.S	Three core Shielded cable	Single conductors	any device
KA1C.1	Three core Shielded cable	Banana plugs	any device
KA1C.2	Three core Shielded cable	7 poles BVT connector	any device
KA1C.3	Three core Shielded cable	Triad01 Palminstruments	Palmsens - Palminstruments
KA1C.4	Three core Shielded cable	BNC connectors	eDAQ
KA1C.5	Three core Shielded cable	Small banana plugs	any device

Recommended Vessel

- TC4, TC5, TC6


Ordering Information

- The order is specified by whole product code
- Minimum order quantity - 1 connector
- Delivery time for standard KA1C connector is 4 weeks from receipt of order
- Delivery time for non-standard KA1C connector depends on final technical specification of order


Example of Order

- 5 pieces - KA1C.S


KA1C.S




KA1C.1




KA1C.2



KA1C.3



KA1C.4



KA1C.5

