

TESTING SENSOR

Type: TS

Description

Electrochemical measurements are very sensitive. Sometimes cables, connectors and their wetness can create the measurement errors. Testing sensor enables simple identification of electrochemical hardware error. The sensor has same geometry as AC1 type but its response produced by resistor. The response is stable and independent on chemistry. Testing sensor is formed on a corundum ceramic base. At the end of the sensor there is a contacting field which is connected with the SMD resistors by the silver conducting paths which are covered by a dielectric protection layer.



Physical Parameters

Dimensions:

Weight: 0.4 gms
 Length: 25.40 mm
 Width: 7.26 mm
 Thickness: 1.2 mm

Evaluation units

- PalmSens
- eDAQ

Connector types for TS sensors range

	KA1	KA1C	KA1s	KA4
TS	✓	✓	✓	✓

Sensor is compatible with FC2, FC4

Sensor Usage

- Simple identification of electrochemical hardware error
- Testing of connector measurement channel for AC1 placing
- Calibration due to measurement of resistor characteristic with defined value of resistance

Ordering information

- The order is specified by whole sensor description formula
- Minimum order quantity - 5 sensors
- All order quantities are to be in multiples of 5 e.g. 10, 15, 20, etc.

Examples of Order

- 10 pieces - TS
- 25 pieces - TS