

## WORKING CLASSIC ELECTRODE CONIC

Type: WCEc.W\*.D\*\*.E\*\*\*

### Description

WCEc is a glass tube with a Pt or Au wire inside. Only circle glossy area of the metal wire is outside the glass tube. Standard cable length of working classic electrode is 1,5 m. Connector or banana plug color is red. Working classic electrode (WCEc) is designed for measurements with TC4, TC5, TC6 glass cell.



### Physical Parameters

#### *Dimensions:*

Length: 57 mm  
 Diameter: 11mm  
 Weight: 50 gms

### Electrode material

WCEc.W\*.D\*\*.E\*\*\* (Working Classical Electrode Conic)

W - Working	3 - platinum wire, $\Phi$ 1000 $\mu$ m
C - Classical	4 - golden wire, $\Phi$ 20 $\mu$ m
E - Electrode	E - Electrode ending
c - Conic	S - Bare wire
W - Electrode material	1 - Banana plug
S - Pt	2 - 7 poles BVT connector
1 - Au	3 - Triad 01 PalmInstruments
D - Diameter of the wire	4 - BNC connector
1 - platinum wire, $\Phi$ 100 $\mu$ m	5 - Small banana plugs
2 - platinum wire, $\Phi$ 600 $\mu$ m	

### Ordering Information


- The order is specified by whole product code
- Minimum order quantity - 1 electrode
- Delivery time for standard WCEc.W\*.D\*.E\* electrode is 4 weeks from receipt of order

- Delivery time for non-standard WCEc.W\*.D\*.E\* electrode depends on final technical specification of order


**Example of Order:** - 10 pieces - WCE.WS.D200.E4

### Internal Wiring

WCEc.W*.D**.ES	WCEc.W*.D**.E1
----------------	----------------



WCEc.W*.D**.E2	WCEc.W*.D**.E3
----------------	----------------



WCEc.W*.D**.E4	WCEc.W*.D**.E5
----------------	----------------

