





Electrochemical tools

Maximize bench space by using to store pipettors when not in use. Holds up to three pipettors. Material of the rack is PMMA.



(Pipettor is not included)

Desiccator

An airtight container which maintains an atmosphere of low humidity through the use of a suitable drying agent which occupies the bottom part of the desiccator. It is for the storage of dry objects that must not be exposed to the moisture normally present in the atmosphere.



Inner Dimensions:

Width 302 mm;

Length 324 mm;

Height 420mm.

Three shelves with 85 mm space.